SDMA VET control kit

For veterinary use only!

Intended use:

For quality control procedures of the quantitative determination with the SDMA VET test kit on a solo or CUBE-VET analyser

 Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

	Order information	Indication	Kit size
	Order number: VT0300 Order number: VT0301 Order number: VT3000	SDMA VET test kit SDMA VET test kit SDMA VET control kit	16 tests 6 tests 1 x 2 ml (decision level)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 – 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

This SDMA VET control kit is a valuable tool for assuring quality and precision of the SDMA determination with the SDMA VET test kit (VT0300 / VT0301) and provides a reliable marker for monitoring the accuracy of results.

Preparation and Handling

The controls are liquid and ready to use. Thoroughly mix the content of the vial before each use by gently inverting it. After use please cap vial immediately and put back into refrigerator. To avoid contamination use single use pipetting tips only.

Stability and Storage

Unopened vials are stable at 2 – 8 °C until the expiration date stated on the label. Once opened, vials are stable for 30 days, when stored tightly capped at 2 – 8 °C. Close immediately after use. Warning: DO NOT FREEZE!

Waste Management

Please refer to local legal requirements

Warnings and Precautions

This control kit is for veterinary diagnostic use only!
DO NOT INGEST! Exercise the normal precautions required for handling all laboratory reagents.

Set the species to „Control“ before running the QC material

Compare the results with the values indicated in the lot-specific value table

If the result is not within the acceptable range, repeat measurement. If result is again not within the acceptable range, please contact your distributor

Durchführung einer SDMA VET Kontrollmessung Processing of a SDMA VET control measurement

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ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) before use!

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2. Control preparation

For further details please see user manual of laboratory photometer

- 2.1 Set species to „Control“
- 2.2 Put a few drops of control liquid into a serum cup by squeezing the bottle.
- 2.3 Aspirate 40 µl control liquid out of serum cup using a pipette

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3. Control processing

- 3.1 Dispense 40 µl control liquid INTO THE LIQUID in the ERS cuvette
- 3.2 Apply ERS cap firmly onto ERS cuvette
- 3.3 Place ERS cartridge into laboratory photometer
- 3.4 Start automatic control processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.

Wertetabelle für das SDMA VET Control Kit Value sheet for the SDMA VET control kit

LOT

2403-1

Kontrollwerte | Control Values

2024-12

Target	value (min)	value (max)
50 µg/dl	34 µg/dl	67 µg/dl
2.48 µmol/l	1.66 µmol/l	3.29 µmol/l

Technical details subject to change without notice.

Doc-ID: 62f6f78f-a1a7-40a9-ad24-0a7b856e0315 ; latest revision: DP 3.0 ; date of print: 2020-09-24

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