T4 VET control kit

For veterinary use only!

For quality control procedures of the quantitative determination with the T4 VET test kit on a solo or CUBE-VET analyser

Eurolyser Diagnostica GmbH
Bindergasse 3
5020 Salzburg/Austria

Tel. +43 (0)662 / 43 21 00
Fax +43 (0)662 / 43 21 00-50
www.eurolyser.com

English

	Order information	Indication	Kit size
	Order number: VT 0140 Order number: VT 0141 Order number: VT 1400	T4 VET test kit T4 VET test kit T4 VET control kit	16 tests/pack 6 tests/pack 1 x 5 ml (decision level)

Test kit preparation: Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

This T4 VET control kit is a valuable tool for assuring quality and precision of the T4 determination with the T4 VET test kit (VT0140/VT0141) and provides a reliable marker for monitoring the accuracy of results.

Required Material

(not supplied with the kit)

- Pipettes

Preparation and Handling

The control is liquid and ready to use. Thoroughly mix the contents of the vial before each use by gently inverting for several minutes. Open the vial and transfer the required quantity of control into a clean sample cup. Use sample cup only for one measurement circle and use new sample cups further on. Replace the cap immediately and store the opened vial at 2 - 8 °C. Always use disposable pipette tips!

Stability and Storage

Unopened vials are stable at 2 - 8 °C until the expiration date stated on the label. Once opened, control vials are stable for 30 days when stored tightly capped at 2 - 8 °C. Bacterial contamination produces an increase in turbidity and/or a characteristic odor. Discard vial if evidence of microbial contamination is observed. **DO NOT FREEZE!**

Waste Management

Please refer to local legal requirements

Warnings and Precautions

This control kit is for veterinary diagnostic use only! **DO NOT INGEST!** Exercise the normal precautions required for handling all laboratory reagents. Components of the control which are derived from human source material have been tested using FDA accepted methods and found nonreactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C (HCV), HIV1 and HIV2. However, no test method can offer complete assurance that products derived from human source material are free of infectious agents.

Set the species to „Control“ before running the QC material

Compare the results with the values indicated in the lot-specific value table

If the result is not within the acceptable range, repeat measurement. If result is again not within the acceptable range, please contact your distributor

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ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20 - 25 °C) before use!

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2. Control preparation

For further details please see user manual of laboratory photometer

- 2.1 Set species to „Control“
- 2.2 Aspirate 20 µl control liquid using a pipette

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3. Control processing

- 3.1 Dispense 20 µl control liquid INTO THE LIQUID in the ERS cuvette
- 3.2 Apply ERS cap firmly onto ERS cuvette
- 3.3 Place ERS cartridge into laboratory photometer
- 3.4 Start automatic control processing by pressing the  start button on the solo laboratory photometer, or by closing the door of the CUBE-VET laboratory photometer.

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Technical details subject to change without notice.

Doc-ID: d68795b0-2f0e-448a-8008-f810ddb79672 ; latest revision: DP 3.0 ; date of print: 2023-10-18

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