














BH073134A








EIGHTCHECK-3WP ASSAY SHEET for pocH-100iV Diff Normal Level



LOT 4111 0822 

Parameters	Mean	Limit %	Limit #	Range	Barcode	
WBC [10 ³ /μL]	7.0	9.0	0.6	7.6		
				6.4		
RBC [10 ⁶ /μL]	4.55	5.0	0.23	4.78		
				4.32		
HGB	[g/dL]	12.0	5.0	0.6		
				12.6		
	[g/L]	120.0	5.0	6.0		
				114.0		
[mmol/L]	7.4	5.0	0.4	7.8		
				7.0		
HCT	[%]	36.3	8.0	2.9		
				39.2		
	[L/L]	0.363	8.0	0.029	0.392	
					0.334	
MCV [fL]	79.8	6.0	4.8	84.6		
				75.0		
MCH	[pg]	26.4	6.0	1.6		
				28.0		
	[amol]	1638	6.0	98	1736	
					1540	
MCHC	[g/dL]	33.1	7.0	2.3		
				35.4		
	[g/L]	331.0	7.0	23.0	354.0	
					308.0	
	[mmol/L]	20.5	7.0	1.4	21.9	
					19.1	

 27-Jul-2024 

Parameters	Mean	Limit %	Limit #	Range	Barcode
PLT [10 ³ /μL]	217	15	33	250	
				184	
LYM% [%]	31.8	13.0	4.1	35.9	
				27.7	
OTHR% [%]	11.5	40.0	4.6	16.1	
				6.9	
EO% ¹⁾ [%]	56.7	10.0	5.7	62.4	
				51.0	
LYM# ¹⁾ [10 ³ /μL]	2.2	20.0	0.4	2.6	
				1.8	
OTHR# ¹⁾ [10 ³ /μL]	0.8	50.0	0.4	1.2	
				0.4	
EO# ¹⁾ [10 ³ /μL]	4.0	15.0	0.6	4.6	
				3.4	

1) These parameters values should be used as a reference only.

